

## ZENEX Fixture Compatibility



Fixture		Abutment	
System	Connection	Company	Connection
I – System	X – Narrow		Yellow 3.0
I – System	Narrow		IS – III Active S – Narrow
			TS Narrow
I – System	Mini		IS Regular
I – System	R/W		IS Regular
T – System	Mini		TS Mini
T – System	R/W		TS Regular
R – System	M/R/W		AnyRidge

## ZENEX Abutment Compatibility



Fixture		Abutment	
Company	Connection	System	Connection
<b>Neo Biotech</b>	IS – III S – Narrow	T – System	Mini
	IS – III Active	I – System	M/R/W
<b>Dentium</b> For Dentists By Dentists	SuperLine II	I – System	M/R/W
<b>OSSTEM<sup>®</sup></b> IMPLANT	TSIII Narrow/Mini	T – System	Mini
	TSIII R	T – System	R/W
<b>MEGAGEN</b>	AnyRidge	R – System	M/R/W
<b>DIO<sup>®</sup> IMPLANT</b>	UFII R/W	I – System	M/R/W
<b>DENTIS</b>	SQ Narrow	I – System	X – Narrow
	SQ R/W	I – System	M/R/W
<b>CWM</b> Cowellmedi Co., Ltd.	INNO Submerged	I – System	M/R/W
	INNO Narrow	T – System	Mini